

TIPS & TECHNIQUES

FOR THE IMPLANTATION OF

ProxiCor[®]
FOR PERICARDIAL CLOSURE

ProxiCor[®]
FOR CARDIAC TISSUE REPAIR

Refer to the Instructions for Use for indications, full implant instructions, warnings and precautions, and contraindications.

DO:

- HYDRATE PROPERLY:** Hydrate (soak) ProxiCor in **room-temperature** sterile saline for **1-2 minutes**. Excessive hydration may result in delamination.
- SIZE MATERIAL APPROPRIATELY:** To promote effective remodeling, use only what is needed of the material. ProxiCor should not be pulled too tight, and sutures should not be pulling. **Minimize material redundancy and tension.**
- USE NON-ABSORBABLE SUTURES:** A non-absorbable, monofilament polypropylene suture is preferred.
- USE EITHER SIDE:** ProxiCor is not sided and you may use either side in your procedure.
- IMPLANT PROPERLY: Suture to healthy, viable tissue** – ProxiCor ECM bioscaffold relies on angiogenesis and cellular infiltration from the sutured edge.
- SUTURE CIRCUMFERENTIALLY:** Make sure all layers of ProxiCor are penetrated to ensure good tissue approximation and regeneration potential.

FOR PROXICOR FOR PERICARDIAL CLOSURE ONLY:

- IMPLANT PROPERLY:** Use a continuous, running stitch (approximately 1 stitch per cm) to approximate the edge.
- FLUSH UNTIL CLEAR:** Irrigate shed blood from pericardial sack with warm saline until clear.
- INSERT CHEST TUBES:** Recommend two drains, one inside and one outside of the reconstructed pericardium.

DON'T:

- USE WITHOUT RECOMMENDED HYDRATION:** Without proper hydration, ProxiCor can tear and not retain sutures.
- USE WITH GLUES, SEALANTS, OR PLATELET GELS:** Avoid using glutaraldehyde-based glue, platelet gels, and hemostatic agents as they may increase crosslinking and negatively impact remodeling.
- SUTURE TO SYNTHETIC OR NON-VIABLE MATERIAL:** Do not sew to synthetic or non-viable material.
- SAVE OR REUSE:** Discard all open or unused portions of ProxiCor.

NOTE:

- ProxiCor is not indicated for the construction or replacement of total valves or conduits.
- Users have reported that when repairing an ASD or VSD, disrupting the endothelial lining where ProxiCor contacts the tissue prior to implantation may facilitate remodeling.
- Use in patients with connective tissue disorders (i.e. Marfan Syndrome) has not been evaluated.

aziyo.com

©2021 Aziyo Biologics, Inc.
All rights reserved.
MK-1467-01/Rev D | US Use Only

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information.

AZIYO[®]
B I O L O G I C S